4^{SENSOR} BSCT

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Operating Instructions

A Rapid Test For β -lactams, Dihydrostreptomycin, Streptomycin, Chloramphenicol and Tetracyclines Detection in MILK

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
1. INTRODUCTION

4Sensor (KIT060) is a rapid qualitative test that allows you to simultaneously detect and discriminate the presence or absence of **β -lactams, Tetracyclines, (Dihydro)-Streptomycin and Chloramphenicol** molecules in a milk sample.

4Sensor is a competitive test involving specific receptors and generic monoclonal antibodies in one single operation. The test requires the use of two components. The first component is a microwell containing predetermined amounts of both receptors and antibodies linked to gold particles. The second is a dipstick made up of a set of membranes with specific capture lines.

When the reagent from the microwell is re-suspended with a milk sample, both receptors and monoclonal antibodies will bind the corresponding analytes if present during the first 5-minute incubation at 40°C. Afterwards, when the dipstick is dipped into the milk, the liquid starts running vertically on the dipstick and passes through capture zones. If the sample is free of antibiotics, a color development occurs at the 'test' lines, indicating the absence of the targeted analytes in the milk sample. On the contrary, the presence of antibiotic(s) in the sample will impair the reagent to color the 'test' line. In any case, the 'control' line must be colored to validate a test.

2. COMPOSITION OF KITS

4Sensor Milk Kits contain everything needed to perform 96 measurements.

- 12 tubes each with 1 strip of 8 reagent microwells and 8 dipsticks.
- 1 micropipette of 200 μ l and disposable tips.
- 1 vial of POSITIVE CONTROL (red cap - the Positive Control is stained by a dye to give it a very pale red color).
- 1 vial of NEGATIVE CONTROL (green cap - the Negative Control is stained by a dye to give it a very pale green color).
- The certificate of compliance.
- The user manual.

Additional equipment/material needed but not included

- 1 HeatSensor (40°C incubation, refer to the used HeatSensor manual).
- 1 ReadSensor (optional, refer to the ReadSensor manual).
- Timer/Watch.
- Distilled water.

3. GENERAL REMARKS

- Store the kit at a temperature between 2°C and 8°C in a controlled fridge.
- Avoid repeated temperature variation.
- Avoid exposure of the reagents to moisture and light.
- Do not mix components of kits from different batches.
- Do not use the kit after the expiration date.
- Avoid touching the membrane on the dipstick.
- Avoid leak of reagent during sample addition and mixing step.
- Empty one tube before opening another tube.

- Close accurately the tube after use.
- The milk sample must be liquid and homogeneous. There can be neither clots or sedimentation phases. The ideal temperature of the milk sample is between 4 and 20°C.
- The temperature to perform the test is 40°C ± 1°C. Use the "HeatSensor". Any other type of incubator is not appropriate to perform the 4Sensor assay. The use of other incubators that indicated in this document is the responsibility of the client and does not bind Unisensor (see « HeatSensor manual »).
- The timing to perform the test is 5 minutes + 5 minutes. This timing must be respected, no extra time (second) is allowed.
- When drying, the color intensities of the lines will become sharper.
- When a positive result is recorded, the test result should be confirmed by performing two more tests.

4. TESTING PROCEDURE

Before Starting

1. Read the instruction carefully.
2. Choose a clean and dry place to perform the test. Wash and dry your hands.
3. Reader and incubator must be clean.
4. Before opening the reagents, take the kit out of the fridge and wait until the temperature of the reagents reaches the room temperature.
5. Switch on the « HeatSensor » and select the program set at the right temperature (40°C) and timing (5 min + 5 min) (see HeatSensor manual). Wait until the temperature has stabilised at 40°C.
6. Switch on the « ReadSensor » if instrumental reading is chosen for result interpretation (see ReadSensor manual); Make sure that the check / calibration was performed according to the User's Manual.

Testing Steps

1. Identify all samples with a number (the positive and negative controls included if necessary).
2. Open a tube and take out as many microwells as there are milk samples to be tested. Be careful to leave the protecting film on the unused microwells. Immediately put the remaining reagents back into the white tube without damaging the dipstick, close the tube and make sure it is tightly sealed.

Testing more than 8 samples at a time is not recommended. In the case of many samples have to be analyzed, try to perform the tests in cascade to avoid any delays occurring during the manipulations. For more than 3 samples, consider using a multichannel pipette. Make sure you have the same incubation time for each sample.

3. For each sample, place a new disposable tip on the micropipette and transfer 200µl of milk into one microwell, then mix gently to homogenize.
4. Place immediately the microwell(s) into the « HeatSensor » (at 40°C), close the device's lid and press the START button on the device. The initial incubation of 5 minutes begins. Place the dipsticks into the HeatSensor holder if using HeatSensor DUO/OCTO.

WARNING: when reagents and milk are in contact, the reaction starts. Make sure the time for sample preparation and handling is minimal.

5. After the initial incubation, the dipstick will be dipped automatically into the microwell if using HeatSensor DUO/OCTO or manually in case of other incubators and incubation will proceed for the next 5 minutes.
6. When the incubation with the dipstick is over, remove the dipstick from the microwells and holder.
7. Gently remove the sample pad without damaging the central membrane.

- Proceed to the interpretation by visual observation or with the ReadSensor (as described below) within 5 minutes.
- If you are not planning to performed any other test in the next 24h, put the tube back for storing into the fridge (2°C to 8°C).

Milk Powder Dilution

In an appropriate bottle, mix 10 g of milk powder with 90 ml of warm (40°C) and distilled water. For an optimal dilution, shake vigorously.

Use Of The Control Vials (OPTIONAL)

- Remove the cap and rubber stopper from the vial and add 1 ml of distilled water.
- Replace the stopper and leave to stand for 1 minute and a half at room temperature.
- Vigorously shake the sample during 10 seconds.
- Check that the sample is properly dissolved.
- The material is now ready to use.

5. INTERPRETATION OF THE RESULTS

Visual Interpretation (See Figure A)

- If the control is not visible (with or without the presence of a test line), do consider the test as not valid, don't interpret the result and directly run another test.
- If the test is valid, examine one test line at a time and compare its color intensity to that of the control line.
 - If the test line is darker in color than the control line, the result is **NEGATIVE**. Given the sensitivity of our test, the milk sample contains no antibiotics or antibiotics at a lower level than the value (limit of detection) mentioned in table A.
 - If the test line's color intensity is similar, lighter (absent) compared to the control line, the result is **POSITIVE**. Given the sensitivity of our test, the milk sample contains antibiotics at or above the value (limit of detection) mentioned in table A.
- When you have interpreted one test line, do the same for the other lines.
- Write down your assessment on each of the dipsticks.

NOTE

- If you hesitate, consider the sample as **POSITIVE** and confirm your interpretation by performing two more tests.
- Dipsticks can be archived as a permanent record if the sample pad has been removed just after the end of the test. Allow the strip to dry before storage. Please note that the color intensity can vary with drying.

ReadSensor Interpretation

1. Make sure you have entered the latest version available of the reading method into your ReadSensor (Check Unisensor Website for the latest appropriate version).
2. If the control is not visible (with or without the presence of test lines), do consider the test as not valid, don't interpret the result and directly run another test. The ReadSensor will confirm it by giving an INVALID result.
3. You should read the dipstick within 5 minutes of performing the test.
4. On the basis of a ratio obtained between the test and control line color's intensities, the reader will interpret the result as follows:

	Results Interpretation	
Test Result :	Negative (NEG)	Positive (POS)
Interpretation :	Contamination absent or below the limit of detection	Contamination at or above the limit of detection

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